

Human Herpes Virus (HHV) Assay Order Form

Please use the following instructions when collecting your samples. It is very important to **read all instructions** before beginning the collection procedure. If you have any questions or problems, do not hesitate to contact us.

Sample collection for HHV real-time PCR assay

Use Universal Precautions when collecting all potentially infectious samples. The required sample is serum, preferably frozen.

1 Collect 7 cc of blood by venipuncture in a serum separator tube. **Do not allow blood to clot prior to centrifugation.**

2 Centrifuge the tube at 3,000 rpm for 10 minutes. The gel should be horizontal in the tube and clearly separating serum and red blood cell fractions. If not, centrifuge an additional 5 minutes.

3 The centrifuged serum tube can be shipped frozen with a cold pack or at ambient temperature.

For questions or concerns contact

Sarah Helber, PhD

Group Leader, Clinical Genetics

Tel 800.735.9224

Fax 804.348.2641

Shipping address

Attn: Clinical Laboratory

AlBioTech

601 Biotech Drive

Richmond, VA 23235

Patient information

Patient name:

Date of birth: / /
 mm / dd / yyyy

Gender: Female Male

Referring physician:

Patient ID#:

Collection date: / /
 mm / dd / yyyy

Volume:

Type of sample (i.e. CSF, serum):

Requested test(s)

Panel (HHV 4, HHV 5, HHV 6, HHV 7)

HHV 1*	HHV 2*	EBV (HHV 4)
CMV (HHV 5)	HHV 6 (a+b)	HHV 6a
HHV 6b	HHV 7	HHV 8

Pricing

	TAT (days)	3-5	5-10
Panel (HHV 4, HHV 5, HHV 6, HHV 7)		\$600	\$500
Individual assays (each)		\$250	\$150

Payment information

Payment must be included with completed order form. Please provide payment in the form of credit card, purchase order or check.

Your sample will not be processed without payment.

*HHV 1 and HHV 2 are eligible for medical reimbursement. Please contact Gynecor® at 866.649.4246.

Name on credit card:

Credit card no.:

Exp. date: /
 mm / yy

Purchase order number:

Check (make check payable to AlBioTech)